


## MICRO Instruments with round Light-weight handle

### *Forceps*




Tip width 0,3mm



**28.0140AL** MICRO-FORCEPS 200mm bay. sharp 0,3mm

Tip width 0,5mm



**28.0142AL** MICRO-FORCEPS 200mm bay. blunt 0,5mm

Tip width 0,2mm



**28.0220AL** MICRO-FORCEPS 15 cm straight sharp

Tip width 0,2mm



**28.0221AL** MICRO-FORCEPS 15 cm curved sharp

Tip width 0,3mm



**28.0230AL** MICRO-FORCEPS 15 cm straight blunt

Tip width 0,3mm



**28.0231AL** MICRO-FORCEPS 15 cm curved blunt

## MICRO Instruments with round Light-weight handle *Forceps*



Tip width 0,3mm

**28.0222AL** MICRO-FORCEPS 18 cm straight sharp



Tip width 0,3mm

**28.0223AL** MICRO-FORCEPS 18 cm curved sharp



Tip width 0,4mm

**28.0232AL** MICRO-FORCEPS 18 cm straight blunt



Tip width 0,4mm

**28.0233AL** MICRO-FORCEPS 18 cm curved blunt

## MICRO Instruments with round Light-weight handle *Forceps*



Tip width 0,3mm



**28.0260AL** MICRO-Tying Forceps 15 cm straight

Tip width 0,3mm



**28.0261AL** MICRO-Tying Forceps 15 cm curved

Tip width 0,6mm



**28.0272AL** MICRO-Tying Forceps 18 cm straight

Tip width 0,6mm



**28.0273AL** MICRO-Tying Forceps 18 cm curved

## MICRO Instruments with round Light-weight handle *Scissors*



Blade length 10 mm



**28.2122AL** MICRO-SCISSORS 15 cm straight

Blade length 10 mm



**28.2123AL** MICRO-SCISSORS 15 cm curved

**28.2123ZAL** MICRO-SCISSORS 15 cm curved  
with one serrated blade

Blade length 10 mm



**28.2122AL** MICRO-SCISSORS 18 cm straight

Blade length 10 mm



**28.2123AL** MICRO-SCISSORS 18 cm curved

## MICRO Instruments with round Light-weight handle *Needleholders*




Tip width 0,8mm



**28.3600AL** MICRO-Needleholder 15 cm straight

Tip width 0,8mm



**28.3601AL** MICRO-Needleholder 15 cm curved

Tip width 0,8mm



**28.3602AL** MICRO-Needleholder 18 cm straight

Tip width 0,8mm



**28.3603AL** MICRO-Needleholder 18 cm curved



## MICRO Instruments with round Light-weight handle ***Needleholders with lock***

Tip width 0,8mm



**28.3630AL** MICRO-Needleholder 15 cm straight w.lock

Tip width 0,8mm



**28.3631AL** MICRO-Needleholder 15 cm curved w.lock

Tip width 0,8mm



**28.3632AL** MICRO-Needleholder 18 cm straight w.lock

Tip width 0,8mm



**28.3633AL** MICRO-Needleholder 18 cm curved w.lock